

Creating add-ins

Data provider, method importer, and integrator contracts; load folders and bootstrap.

Quanto loads **three kinds of add-ins** at startup from folders next to the app. Drop a DLL in the right folder; no registry keys.

Kind	Folder	Contract	Job
Data provider	AddIns\DataProvider\	IDataProvider	Find samples, read bottle metadata, extract chromatograms
Method importer	AddIns\MethodImporter\	IMethodImporter	Find and parse vendor methods into compounds / signals / levels
Integrator	AddIns\Integrator\	IIntegrator	Detect peaks on <code>x[]</code> / <code>y[]</code>

Contracts live in `Quanto.Shared.Core` (`Quanto.Shared.Core.Interfaces`). The host loader is `Quanto.AddInFactory.AddInLoader` .

How loading works

1. At startup the host calls `AddInLoader.Load(integratorPath, dataProviderPath, methodImporterPath)` .
2. Each folder is scanned for `*.dll` . Assemblies are `Assembly.LoadFrom` .
3. Instantiation is **lazy** — constructors (and heavy vendor init) run on first use, not at discovery.
4. Registration key is the **concrete class name** (e.g. `AgilentMassHunterDataProvider` , `Quanto` , `SciexWiffMethod`).
5. Optional `IAddIn: AlternativeDllPath` (vendor DLL search folder) and `Version` (for logs).
6. Optional **bootstrap** class: a tiny type with an `AlternativeDllPath` property (conventionally named `AddInBootstrap` , or any `*Bootstrap`). Instantiated first so vendor paths are registered **before** reflecting types that reference vendor assemblies.

Paths (from `AppPaths`):

```
text {app}\AddIns\ DataProvider\ MethodImporter\ Integrator\ VendorSupportFiles\Agilent\
VendorSupportFiles\Sciex\
```

Relative `AlternativeDllPath` values are resolved from the application directory (examples from shipped add-ins: `AddIns\VendorSupportFiles\Agilent` , `AddIns\VendorSupportFiles\Sciex`).

Shipped projects copy only their **own** primary DLL into the staging add-in folder; vendor natives stay under `VendorSupportFiles\...` .

Project setup

- Target **Windows x64** and the same Quanto Windows TFM as the host (.NET 10 Windows).
- Reference `Quanto.Shared.Core` only for contracts and shared types/enums (do not take a hard dependency on `Quanto.WPF`).
- Implement the matching interface(s) on a **public non-abstract class**.

- Optionally implement `IAddIn`.
- Add a dependency-free bootstrap type if you need vendor DLLs on the probe path before your main type loads.
- Build output: one add-in DLL into the correct `AddIns\...` subfolder.

Example bootstrap (duck-typed; implementing `IAddInBootstrap` is fine but not required for discovery):

```
```csharp
namespace MyVendorData;

public sealed class AddInBootstrap {
 public string? AlternativeDllPath =>
 @"AddIns\VendorSupportFiles\MyVendor";
}
```

## Data provider

`IDataProvider` extends:

- `ISampleFinder` — `HashSet<SampleFile> FindSamples(string path)`
- `ISampleFileReader` — `void ReadSampleFile(SampleFile sampleFile)`
- `ISignalReader` — chromatogram extraction

```
```csharp
public interface ISignalReader {
    SignalPoints GetSignalPoints(
        string filePath, double rt, double
        beginOffsetTime, double endOffsetTime, double mass1, double beginOffsetMass1, double
        endOffsetMass1, double mass2, double beginOffsetMass2, double endOffsetMass2, IonPolarity?
        ionPolarity = null, double? fragmentorVoltage = null, double? collisionEnergy = null,
        ChromatogramAcquisitionMode preferredMode = ChromatogramAcquisitionMode.Auto);
}
```

```
SignalPoints GetTotalSignal(string filePath);
// optional override:
SignalPoints GetTotalSignal(string filePath, ChromatogramAcquisitionMode mode);
}
```

`ChromatogramAcquisitionMode: Auto, Scan, Sim, MsMs.`

`FindSamples` should return `SampleFile` instances with `DataProvider` set to this instance and fill what you can (name, AcqDateTime, instrument, position, dilution, dose hints). The host matches sample type / dose from the method afterward.

Optional host hooks (implement if useful): `IMrmTransitionCatalog`, `IDataFilePrefetch`, `IDataFileRelease`.

Reference implementations

Class	Project	Format
<code>AgilentMassHunterDataProvider</code>	<code>Quanto.AddIn.DataProvider.AgilentMHDData</code>	MassHunter .D
<code>AgilentChemStationDataProvider</code>	<code>Quanto.AddIn.DataProvider.AgilentCSData</code>	ChemStation .D (native MSDChem)
<code>SciexWiff</code>	<code>Quanto.AddIn.DataProvider.SciexWiff</code>	SCIEX .wiff

Method importer

```
csharp public interface IMethodImporter { string Name { get; set; } HashSet<MethodFile> FindMethods(string path); void LoadMethod(string filePath); void ReadMethodFile(MethodFile info); ICompound[] GetCompounds(); }
```

Map vendor content onto `ICompound` / `ISignal` / `ILevel`:

- Compound: name, group, expected time, begin/end offsets, ISTD flag and link, curve fit/weight/origin, concentration limits, signals, levels.
- Signal: name, `Mass1` / `Mass2`, collision energy, fragmentor, ion polarity, `AcquisitionMode`, quant/qualifier flags, % and tolerance.
- Optional `ICompoundWithAdditionalQuants` when the vendor has extra quant slices (e.g. MassHunter RA variants).

`FindMethods` should attach `MethodImporter = this` on each `MethodFile` and set a display `Name`.

Reference implementations

Class	Project	Source
<code>AgilentMhMethod</code>	<code>Quanto.AddIn.MethodImporter.AgilentMHMethod</code>	<code>.m</code> / <code>.quantmethod.xml</code>
<code>AgilentMassHunterDataMethod</code>	<code>Quanto.AddIn.MethodImporter.AgilentMHData</code>	Method from <code>.D</code>
<code>AgilentCsMethod</code>	<code>Quanto.AddIn.MethodImporter.AgilentCSMethod</code>	<code>.M</code> + <code>qdb.mth</code>
<code>SciexWiffMethod</code>	<code>Quanto.AddIn.MethodImporter.SciexWiff</code>	Method from <code>.wiff</code>

Integrator

```
csharp public interface IIntegrator { PeakValues[] Integrate(double[] x, double[] y, string settings); }
```

`PeakValues` fields: `BegX`, `EndX`, `BegB`, `EndB` (peak time bounds and baseline height at those bounds).

`settings` is the free-form string from the channel's integrator settings. Parse what you need; ignore unknown keys.

Channels pick an integrator **by class name** (default shipped name: `Quanto`).

What the integrator returns vs what Quanto calculates

An integrator **only detects peak boundaries**. It does **not** compute area, height, apex time, width, or symmetry.

After your add-in returns `PeakValues[]`, the host always runs **Quanto's shared peak-metric methods** in `Quanto.Shared.Core.PeakCalculations` (`IPeakMetricsCalculator` / `PeakMetricsCalculator`) on the same `x[]` / `y[]` and your `BegX` / `EndX` / `BegB` / `EndB`. That path is used for:

- automatic integration (`ImportPeakValuesService` after any integrator, including custom ones)
- manual re-integration and boundary edits (`PeakIntegrationService`)

Derived metrics include (see `PeakMetrics`):

Metric	Meaning
Area	Trapezoidal integral of the baseline-corrected signal
Height	Apex height above the linear baseline
CenterX / Time	Apex retention time
MaxSignal	Apex signal value
FwAt50 , FwAt5	Peak widths at 50% and 5% of height
Symmetry	Peak symmetry from those crossings
BegR , EndR	Residuals / edge refinements from the shared calculator

Response for quantitation is then **area** or **height** according to the analyte's `QuantitateByHeight` flag — still from these shared metrics, not from the integrator DLL.

So when you write an integrator:

- 1. Focus on finding good start/end times and baseline points.
- 2. Return accurate `BegX` , `EndX` , `BegB` , `EndB` .
- 3. **Do not** reimplement area/height in the add-in; Quanto will calculate them consistently for every integrator (Quanto, Slice, or yours).
- 4. You may call `PeakMetricsCalculator.Instance` yourself for diagnostics or unit tests (same algorithms), but the host will recalculate when importing peaks into the batch.

Reference implementations

Class	Project	Notes
Quanto	<code>Quanto.AddIn.Integrator.Quanto</code>	Production peak detector (boundaries only; metrics still via shared calculator)
Slice	<code>Quanto.AddIn.Integrator.Slice</code>	Minimal stub (returns no peaks) — useful as a template

Minimal skeletons

Integrator

```
```csharp using Quanto.Shared.Core.Interfaces; using Quanto.Shared.Core.Types; // Optional for offline checks: // using Quanto.Shared.Core.PeakCalculations; namespace MyIntegrator; public sealed class MyIntegrator : IIntegrator, IAddIn { public string? AlternativeDllPath => null; public string? Version => "1.0.0"; public PeakValues[] Integrate(double[] x, double[] y, string settings) { // Detect peaks; return only BegX/EndX/BegB/EndB. }
```

```
// Host applies PeakMetricsCalculator for Area, Height, Fw*, Symmetry, etc.
return [];
}
}'''
```

## Data provider / method importer

Follow the same pattern: implement the interface + optional `IAddIn`, discover files under the path the host passes, and return shared types from `Quanto.Shared.Core.Types / Enums`. Prefer looking at the matching Agilent or SCIEX project above rather than re-deriving vendor binary formats.

## Checklist before shipping an add-in

- [ ] DLL builds x64 / Quanto Windows TFM and copies to the correct `AddIns\...` folder
- [ ] Class is public, non-abstract, and assignable to the right interface
- [ ] Class **name** is unique among loaded add-ins of that kind (registry key)
- [ ] Bootstrap + `AlternativeDllPath` set if vendor natives are required
- [ ] Vendor DLLs installed under `VendorSupportFiles\...` (not duplicated beside every add-in)
- [ ] Constructors avoid crashing the process at first use (lazy load still runs your ctor once)
- [ ] Version visible via `IAddIn.Version` or assembly informational version
- [ ] Smoke-test: host log shows Discovered / instantiated lines; Find / Load / Integrate works on a real file
- [ ] Integrators: confirm Area/Height appear after Analyze (shared `PeakMetricsCalculator`), not only raw boundaries

## Related

- Formats and modes: [Compatibility](#)
- Install tree: [Installation and folders](#)
- Analyst import UI: [Reports, import, and export](#)