

Silent deploy

/VERYSILENT checklist, .NET 10 Desktop x64, no post-install launch, TDA optional.

Checklist for SCCM / Intune / unattended install from the Inno Setup section of [Installation and folders](#).

Prerequisites

Item	Notes
OS	Windows 10/11 64-bit
Runtime	.NET 10 Desktop Runtime (x64) unless your package already bundles it
Elevation	Run the setup elevated
TDA converter	Optional — only for MassHunter .D that lack MSTree2.bin (Agilent MassHunter Quant). ChemStation and SCIEX do not need it

Typical silent command

Example setup name (match your build artifact):

```
bat Quanto-1_8-setup.exe /VERYSILENT /SUPPRESSMSGBOXES /NORESTART /LOG="%TEMP%\quanto-setup.log"
```

Parameter	Purpose
/VERYSILENT	No prompts, no progress window
/SUPPRESSMSGBOXES	Auto-answer message boxes
/NORESTART	Do not reboot even if Windows asks
/LOG=...	Detailed Inno log

Also available: /SILENT (progress window), /DIR= , /GROUP= , /CLOSEAPPLICATIONS , /FORCECLOSEAPPLICATIONS , /NOCLOSEAPPLICATIONS . There are no optional Tasks/Components in this script.

After silent install

- The setup does **not** launch Quanto when /SILENT or /VERYSILENT was used (skipifsilent). Start the app yourself or via a shortcut / software-center step.
- Install location is C:\Program Files\Quanto unless /DIR= overrides it.
- User profile folders are created on first run under %LOCALAPPDATA%\Quanto .

Permissions pointer

Use the permissions table in [Installation and folders](#): Program Files read/execute; Public Documents license/config as needed; Local AppData full control; study folders and Agilent .D write when MH indexing must run on the PC.

Related

- [Licensing](#) — shared license folder after deploy
- [Validation notes](#)
- [Troubleshooting](#) — add-ins / TDA